
## **InPro® 6800 Series O<sub>2</sub> Sensors**

**Instruction manual  
Bedienungsanleitung  
Instructions d'utilisation**

**English** **Page** **3**

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# **InPro<sup>®</sup> 6800 Series O<sub>2</sub> Sensors**

## **Instruction manual**

## Contents

|           |   |           |
|-----------|---|-----------|
| <b>1</b>  | <b>Introduction</b> .....   | <b>5</b>  |
| <b>2</b>  | <b>Important notes</b> .....  | <b>6</b>  |
| 2.1       | Notes on operating instructions .....   | 6         |
| 2.2       | Intended use .....  | 6         |
| 2.3       | Safety instructions .....   | 7         |
| 2.4       | Examples of some typical applications .....   | 8         |
| 2.5       | Use in Ex-zones.....  | 8         |
| 2.6       | Ex-classification ATEX .....  | 9         |
| 2.6.1     | Introduction.....   | 9         |
| 2.6.2     | Rated data .....  | 9         |
| 2.6.3     | Special conditions.....   | 10        |
| 2.7       | Ex-classification FM approved .....   | 11        |
| 3.1       | General information.....  | 12        |
| 3.2       | Principle .....   | 12        |
| 3.3       | Scope of delivery .....   | 12        |
| 3.4       | Equipment features .....  | 13        |
| <b>4</b>  | <b>Installation</b> .....   | <b>15</b> |
| 4.1       | Mounting the sensor.....  | 15        |
| 4.2       | Connection.....   | 15        |
| 4.2.1     | Connecting the InPro 6800 to a VP cable .....   | 15        |
| 4.2.2     | Connecting the VP cable to the transmitter .....  | 16        |
| <b>5</b>  | <b>Operation</b> .....  | <b>17</b> |
| 5.1       | Start-up and polarizing .....   | 17        |
| 5.2       | Calibration .....   | 18        |
| 5.2.1     | Purpose of calibration .....  | 18        |
| 5.2.2     | What you have to know for calibration .....   | 18        |
| 5.2.3     | Single point calibration .....  | 19        |
| 5.2.4     | Dual point calibration .....  | 20        |
| <b>6</b>  | <b>Maintenance</b> .....  | <b>21</b> |
| 6.1       | Inspection of the sensor.....   | 21        |
| 6.1.1     | Visual inspection .....   | 21        |
| 6.1.2     | Testing the sensor with the METTLER TOLEDO<br>O <sub>2</sub> Sensor-Master InPro 6800 ..... | 22        |
| 6.1.3     | Testing the sensor via a transmitter .....  | 23        |
| 6.2       | Changing the electrolyte, the membrane body or<br>the interior body.....                    | 24        |
| <b>7</b>  | <b>Storage</b> .....  | <b>27</b> |
| <b>8</b>  | <b>Product specification</b> .....  | <b>27</b> |
| 8.1       | Certificates .....  | 27        |
| 8.2       | Specifications .....  | 28        |
| <b>9</b>  | <b>Ordering information</b> .....   | <b>29</b> |
| 9.1       | Sensors .....   | 29        |
| 9.2       | Accessories.....  | 29        |
| 9.3       | Spare parts .....   | 30        |
| 9.4       | Recommended transmitters.....   | 30        |
| 9.5       | Recommended housings .....  | 30        |
| <b>10</b> | <b>Theory of the polarographic sensor</b> .....   | <b>31</b> |
| 10.1      | Introduction.....   | 31        |
| 10.2      | Principle of the design of an oxygen electrode .....  | 32        |
| 10.3      | Parameters determining current .....  | 32        |
| 10.4      | Polarization voltage .....  | 33        |
| 10.5      | Temperature .....   | 33        |
| 10.6      | Dependence on flow.....   | 34        |
| 10.7      | Oxygen partial pressure – oxygen concentration .....  | 35        |

## 1 Introduction

Thank you for buying the **InPro 6800 sensor** from **METTLER TOLEDO**.

The construction of the InPro 6800 sensors employs leading edge technology and complies with safety regulations currently in force. Notwithstanding this, improper use could lead to hazards for the user or a third-party, and/or adverse effects on the plant or other equipment. **Therefore, the operating instructions must be read and understood by the persons involved before work is started with the sensor.**

The instruction manual must always be stored close at hand, in a place accessible to all people working with the InPro 6800.

If you have questions, which are not or insufficiently answered in this instruction manual, please contact your METTLER TOLEDO supplier. They will be glad to assist you.

## 2 Important notes

### 2.1 Notes on operating instructions

These operating instructions contain all the information needed for safe and proper use of the InPro 6800 sensor.

The operating instructions are intended for personnel entrusted with the operation and maintenance of the sensors. It is assumed that these persons are familiar with the equipment in which the sensor is installed.

#### Warning notices and symbols

This instruction manual identifies safety instructions and additional information by means of the following symbols:

This symbol draws attention to **safety instructions and warnings of potential danger** which, if neglected, could result in injury to persons and/or damage to property.

This symbol identifies **additional information and instructions** which, if neglected, could lead to defects, inefficient operation and possible loss of production.

### 2.2 Intended use

Mettler-Toledo InPro 6800 sensors are intended solely for inline measurement of the oxygen partial pressure in liquids and gases, as described in this instruction manual.

Any use of these sensors which differs from or exceeds the scope of use described in this instruction manual will be regarded as inappropriate and incompatible with the intended purpose.

The manufacturer/supplier accepts no responsibility whatsoever for any damage resulting from such improper use. The risk is borne entirely by the user/operator.

Other prerequisites for appropriate use include:

- compliance with the instructions, notes and requirements set out in this instruction manual.
- acceptance of responsibility for regular inspection, maintenance and functional testing of all associated components, also including compliance with local operational and plant safety regulations.
- compliance with all information and warnings given in the documentation relating to the products used in conjunction with the sensor ( housings, transmitters, etc.).
- observance of all safety regulations governing the equipment in which the sensor is installed.

- correct equipment operation in conformance with the prescribed environmental and operational conditions, and admissible installation positions.
- consultation with Mettler-Toledo Process Analytics in the event of any uncertainties.

## 2.3 Safety instructions

- The plant operator must be fully aware of the potential risks and hazards attached to operation of the particular process or plant. The operator is responsible for correct training of the workforce, for signs and markings indicating sources of possible danger, and for the selection of appropriate, state-of-the-art instrumentation.
- It is essential that personnel involved in the commissioning, operation or maintenance of these sensors or of any of the associated equipment (e.g. housings, transmitters, etc.) be properly trained in the process itself, as well as in the use and handling of the associated equipment. This includes having read and understood this instruction manual.
- The safety of personnel as well as of the plant itself is ultimately the responsibility of the plant operator. This applies in particular in the case of plants operating in hazardous zones.
- The oxygen sensors and associated components have no effect on the process itself and cannot influence it in the sense of any form of control system.
- Maintenance and service intervals and schedules depend on the application conditions, composition of the sample media, plant equipment and significance of the safety control features of the measuring system. Processes vary considerably, so that schedules, where such are specified, can only be regarded as tentative and must in any case be individually established and verified by the plant operator.
- Where specific safeguards such as locks, labels, or redundant measuring systems are necessary, these must be provided by the plant operator.
- A defective sensor must neither be installed nor put into service.
- Only maintenance work described in this operating instruction may be performed on the sensors.
- When changing faulty components, use only original spare parts obtainable from your METTLER TOLEDO supplier (see spare parts list, «Section 9.3»).

- No modifications to the sensors and the accessories are allowed. The manufacturer accepts no responsibility for damages caused by unauthorised modifications. The risk is borne entirely by the user.

## 2.4 Examples of some typical applications

Below is a list of examples of typical fields of application for the oxygen sensors. This list is not exhaustive.

### Measurement in liquids:

- Fermentation
- Yeast propagation
- Wort aeration
- Spring water conditioning
- Storage and processing of fruit juices

### Measurement in gases:

- Monitoring of the oxygen limit value to protect products from oxidation.
- Monitoring of the oxygen limit value in the exhaust air of fermentation vessels.
- Monitoring of oxygen limit values during inertization processes.
- Measurement of oxygen concentrations in process gas mixtures.

## 2.5 Use in Ex-zones

### Attention!

For an installation in Ex-zones please read the guidelines following hereafter:

### Ex-classification ATEX:

II 1/2G Ex ia IIC T6/T5/T4/T3

II 1/2D Ex iaD A20/21

T 69 °C/T 81 °C/T 109 °C/T 161 °C

### Marking and number of the test certificate:

SNCH 01 ATEX 3277 X

### Ex-classification FM approved:

IS/I, II, III/1/ABCDEFGG/T6 Ta = 60 °C

- 53 800 002; Entity

## 2.6 Ex-classification ATEX

### 2.6.1 Introduction

According to RL 94/9/EG (ATEX 95) Appendix I, InPro 6XXX\*/\*\*/\*/\* oxygen sensors are devices group II, category 1/2G and according to RL 99/92/EG (ATEX 137) may be used in zones 0/1 or 0/2 and gas groups IIA, IIB and IIC that are potentially explosive due to combustible substances in the temperatures classes T3 to T6.

For use/installation, the requirements of EN 60079-14 must be observed.

According to RL 94/9/EG (ATEX 95) Appendix I, InPro 6XXX\*/\*\*/\*/\* oxygen sensors are devices group II, category 1/2D and according to RL 99/92/EG (ATEX 137) may also be used in zones 20/21 that contain combustible dusts.

For use/installation, the requirements of EN 61241-14 must be observed.

The sensor measuring circuit, the temperature measuring circuit and the data chip circuit are part of a common intrinsically safe system and are for operation connected to a separately certified transmitter.

The sensor measuring circuit, the temperature measuring circuit and the data chip circuit as part of an intrinsically safe system are isolated from the not-intrinsically safe electric circuits up to a maximum rated voltage of 375 V and from grounded parts up to a maximum rated voltage of 30 V.

### 2.6.2 Rated data

#### Measuring circuit:

Fail-safe ignition protection class Ex ia IIC only for connection to a certified fail-safe circuit.

#### Maximum values:

$$U_i \leq 16 \text{ V}, I_i \leq 190 \text{ mA}, P_i \leq 200 \text{ mW}$$

$$L_i = 0 \text{ (effective internal inductivity)}$$

$$C_i = 900 \text{ pF (effective internal capacitance)}$$

#### Note:

The above maximum values are each the total of all individual circuits of the associated intrinsically safe power supply and transmitter.

### 2.6.3 Special conditions

- The maximum permissible environment resp. medium temperatures for the zone 0 (combustible gases or combustible liquids) are in accordance with the temperature classes shown in the table below:

| <b>Temperature class</b> | <b>Max. environment resp. media temperature</b> |
|--------------------------|---|
| T 6                      | 68 °C   |
| T 5                      | 80 °C   |
| T 4                      | 108 °C  |
| T 3                      | 160 °C  |

- The maximum surface temperatures for the zone 20 (combustible dusts) are in accordance with the environment resp. medium temperatures shown in the table below:

| <b>Surface temperature</b> | <b>Max. environment resp. media temperature</b> |
|----------------------------|---|
| T 69 °C                    | 68 °C   |
| T 81 °C                    | 80 °C   |
| T 109 °C                   | 108 °C  |
| T 161 °C                   | 160 °C  |


Additionally the requirements from section 6 of EN 50281-1-2 have to be followed.

- The capacitance and inductance of the connecting cable must be taken into account in the design.
- The oxygen sensors (O<sub>2</sub> sensors) can be used in/with the housings InFit 76\*-\*-\* resp. InTrac 7\*\*-\*-\* or in/with other suitable housings in hazardous areas.
- The metal body of the O<sub>2</sub> sensors resp. the safety weld-in-sockets resp. the independent housing are, if necessary, to be included into the periodic pressure test of the unit.
- The metal body of the O<sub>2</sub> sensors resp. the safety weld-in-sockets resp. the independent housing must be electrically connected to the potential equalizing system of the installation.

2.7

Ex-classification FM approved

## 3 Product description

### 3.1 General information

The **sensor InPro 6800** with integrated temperature probe is used for oxygen measurement.

The sensor is **sterilizable and autoclavable** and **compatible with CIP** (cleaning in place).

InPro 6800 sensors with ISM functionality offer Plug and Measure as well as enhanced diagnostics features.

ISM sensors are available with K8S (fully digital) or VP connector.

### 3.2 Principle

Here is a short summary of the principle of polarographic measurement on which this sensor is based (Clark 1961).

- a) The Clark polarographic sensor basically consist of a working electrode (cathode), a counter/reference electrode (anode), and an oxygen permeable membrane which separates the electrodes from the sample medium.
- b) The transmitter supplies a constant polarization voltage to the cathode, needed to reduce oxygen.
- c) The oxygen molecules which migrate through the permeable membrane are reduced at the cathode. At the same time, oxidation takes place at the anode and oxidized anode metal (silver) is liberated as silver ions into the electrolyte. The electrolyte closes the electric circuit between anode and cathode (ion conductivity).
- d) The current produced by the reactions described above is measured by the transmitter and is proportional to the partial pressure of oxygen ( $pO_2$ ) in the sample medium.

Please refer to Section 10 «Theory of the polarographic sensor» for further information.


### 3.3 Scope of delivery


Each sensor is supplied fully assembled and factory-tested for correct function together with:

- an electrolyte bottle (34 100 2016)
- a quality control certificate
- inspection certificates 3.1 B (complying with EN 10204.3/1 B)

### 3.4 Equipment features

#### 12 mm Sensor

**25 mm Sensor**

## 4 Installation

### 4.1 Mounting the sensor

**Important!** Remove the green protective cap before mounting the sensor.

#### Mounting the sensor in a housing

Please refer to the instruction manual of your housing explaining on how to mount the sensor in place.

#### Mounting the sensor directly on a pipe or a vessel

The **12 mm sensors** can be mounted directly through a socket with inside thread Pg 13.5 and securely tightened via the Pg 13.5 threaded sleeve.

The **25 mm sensors** can be mounted directly through a standard weld-in socket or the safety weld-in socket from METTLER TOLEDO.

### 4.2 Connection


#### 4.2.1 Connecting the InPro 6800 to a cable


Sensors with ISM functionality require the use of a special AK9 cable as well as an ISM-compatible O<sub>2</sub> transmitter. The sensor is connected to the transmitter via a AK9 cable. The AK9 cable ensures a secure connection between the transmitter and the sensor under harsh industrial conditions. The robust watertight IP68 connector housing guarantees maximum process safety. If you have chosen to use a sensor with digital ISM functionality, the K8S connector will incorporate a pre-amplifier. It is absolutely necessary to protect this electronic component against any electrical charge.

**Do not touch the sensor at the VP connector plug.** To connect the VP cable to the sensor align the slit of the VP connector with the pin in the plug. Then tightly screw the plug to fasten the two parts.

#### 4.2.2 Connecting the VP- or AK9 cable to the transmitter

 **Note:** Cable assignment can be found in the METTLER TOLEDO VP or AK9 cable instruction manual.

 **Note:** For connecting the cable to the terminals of the transmitter, please refer to the instructions given in the METTLER TOLEDO transmitter manual.

## 5 Operation

**Attention!** Before using the sensors for the first time, the electrolyte should be replaced (see «Chapter 6.2»). Due to possible adverse conditions during transport and storage (e.g. airfreight; pressure and temperature variations), the quality of the electrolyte may become impaired. Poor electrolyte quality can lead to erroneous measurement values.

### 5.1 Start-up and polarizing

**Attention!** The protective cap at the tip of the sensor must be removed before putting the sensor into operation.

When the system is operated for the first time or if the sensor has been disconnected from the voltage source (transmitter or O<sub>2</sub> Sensor-Master) for longer than 5 minutes, the sensor has to be polarized prior to calibration by connecting it to the operating O<sub>2</sub> transmitter or to a sensor master. After 6 hours, the sensor is fully polarized and ready for operation.

A shorter polarization period is sufficient if the sensor has been disconnected for only a few minutes. The following table serves to establish the correct polarization time in relation to the depolarization time.

| Depolarization time <sup>1)</sup><br>t <sub>depol</sub> [Min.] | Minimum required<br>polarization time <sup>2)</sup> [Min.] |
|--|--|
| t <sub>depol</sub> > 30  | 360  |
| 30 > t <sub>depol</sub> > 15                                   | 6 * t <sub>depol</sub>                                     |
| 15 > t <sub>depol</sub> > 5                                    | 4 * t <sub>depol</sub>                                     |
| t <sub>depol</sub> < 5   | 2 * t <sub>depol</sub>                                     |

- 1) Depolarization time:  
Time span in which the polarization voltage is cut off from the sensor. This is the case during:
- change of electrolyte
  - change of membrane body
  - the time the cable is disconnected or no transmitter or sensor master is connected to the cable
- 2) Polarization time:  
Time span during which the sensor is under a polarization voltage.

**Attention!** Setting of the polarization voltage on the transmitter for correct measurements:

- Standard applications (e.g. measurement in biotechnology): **-675 mV**
- Measurement of permanently low oxygen concentrations <500 ppb in the presence of volatile acidic components (e.g. carbon dioxide during measurements in breweries): **-500 mV**

**Important!** To ensure the supply of the correct polarization voltage the transmitter must be set accordingly.

## 5.2 Calibration

### 5.2.1 Purpose of calibration

Each oxygen sensor has its own individual slope and own individual zero point. Both values are subject to change, for example, through electrolyte consumption or after exchange of electrolyte or membrane body. To ensure high measurement accuracy of the sensor, a calibration must be carried out after each change of electrolyte or membrane. Prior to calibration, the sensor has to be polarized for at least 6 hours.

**A zero point calibration is only advisable if very high accuracy is required at low oxygen concentrations.**


**To check if your sensor needs a recalibration, you may dry it and take it in the air to check that the reading is close to 100%. If not, then the sensor needs a new calibration.**

### 5.2.2 What you have to know for calibration

**Attention!** If dissolved oxygen is being measured at low concentrations (<500 ppb) in the presence of volatile acidic components (e.g. carbon dioxide (CO<sub>2</sub>) measurements in breweries), the calibration procedure described below should be followed in order to achieve the best measurement performance.

No special calibration procedure is required if dissolved oxygen is being measured in standard applications (e.g. measurement in biotechnology).

**First step:** Adjust the transmitter polarization voltage to  $-675$  mV. Detailed information is given in the transmitter instruction manual.

**Second step:** After adjusting the polarization voltage to  $-675$  mV, equilibration time is required. Wait for about 5 minutes before carrying out step 3.

**Third step:** Perform calibration according to the transmitter' instruction manual.

**Fourth step:** Adjust the transmitter polarization voltage back to  $-500$  mV.

**Steps 1 to 4** can be accomplished automatically if you use the M 700 transmitter. Thus the routine can be minimized to a few keystrokes. In addition the M 700 transmitter must be equipped with the software function SW-700-011 «high CO<sub>2</sub> compensation». Please ask your local Mettler-Toledo distributor.

#### **General remarks:**

- **For calibration in air, the sensor membrane must be dry**, since adhering water drops can falsify the measured oxygen value.
- Make sure that the **oxygen saturation index** of the calibration medium is **correct and remains constant** during calibration.
- In the event of calibration in water or sample medium, the **calibration medium must be in equilibrium with the air**. Oxygen exchange between water and air is only very slow. Therefore it takes quite long time until water is saturated with atmospheric oxygen.
- For correct calibration, a minimum flow rate of the calibration medium is necessary.
- **Calibration in a fermenter should be performed after sterilization** (as sterilization may alter the sensor slope), but prior to inoculation. If it is not possible to perform the calibration after sterilization, the use of an existing membrane body that has been pre-sterilized while mounted in the sensor is recommended. A slope alteration of some per cent can occur with new membrane bodies, particularly after a first sterilization, as the tension of the membrane may be altered by the sterilization process.
- Make sure that all other parameters, such as temperature and pressure, are constant.

For continuous applications, we recommend **periodic recalibration in line with your requirements on accuracy, the type of process in operation and your own experience**. The frequency of the need for recalibration depends very much on the specific application, and therefore appropriate intervals cannot be exactly defined here.

### **5.2.3 Single point calibration**

By carrying out a single point calibration, the factual slope of the sensor can be established. The calibration medium can be water with known oxygen saturation index (e.g. air-saturated water) or air with known water-vapor saturation (e.g. water-vapor saturated air).

After the sensor signal has stabilized, the complete measuring system can then be calibrated to the 100 % value of the desired measurable variable, e.g. 100 %

air, 20.95 % O<sub>2</sub>, or 8.26 ppm at 25 °C (77 °F) and normal pressure (see instruction manual for the transmitter).

Single point calibration should be sufficient for almost all process applications.

#### 5.2.4 Dual point calibration

By carrying out a dual point calibration both slope and zero point of the sensor can be established.

**Important!** In case of a **dual point calibration, always start by the zero point calibration** before calibrating the slope.

Due to the very low zero current of METTLER TOLEDO sensors, dual point calibration is normally not necessary for standard applications.

As a rule, the zero point should be adjusted to zero manually, or it is automatically performed by the transmitter (see instruction manual for the transmitter). A zero point calibration is only advisable if very high accuracy is required at low oxygen concentrations.

**Attention! Incorrect zero point calibration is a frequent source of measurement error. For correct calibration, we recommend the use of nitrogen gas or other oxygen-free medium with a level of purity of at least 99.995 %.**

After the sensor signal has stabilized (after 20...30 minutes), the sensor can be calibrated through the relevant transmitter to the 100 % value of the desired measurable variable, e.g. 0 % air, 0.0 % O<sub>2</sub>, or 0.0 ppm (see instruction manual for the transmitter).

## 6 Maintenance

### 6.1 Inspection of the sensor

#### 6.1.1 Visual inspection

To check your sensor, we recommend the following procedure:

- The contacts of the connector must be dry. Moisture, corrosion and dirt in the connector can lead to false readings.
- Check the cable for buckling, brittle areas or ruptures.
- Before calibration always examine the membrane foil optically for signs of damage. The foil must be intact and clean. Dirty membranes should be wiped clean using a soft, moist tissue.  
Note: An undulated membrane has no influence on the sensor performance, assuming the membrane is intact.
- The membrane body must be replaced if the sensor has too long a response time, the reading is unstable or subject to drift, and if the sensor cannot be calibrated or the membrane shows sign of mechanical damage.
- Check the cathode area for discoloration, contamination or cracks in the glass. If necessary rinse with demineralized water and clean with a clean soft brush or soft paper tissue.

**Attention! Do not use any cleaning agents containing alcohol. This could damage the sensor or lead to fault current.**


**Attention! The glass body is fragile and sensitive to vibration.**

### 6.1.2 Testing the sensor with the METTLER TOLEDO O<sub>2</sub> Sensor-Master InPro 6800

We recommend to use the METTLER TOLEDO O<sub>2</sub> Sensor-Master InPro 6800 to check the quality of your sensor as follows (not applicable for digital ISM sensors):

- Connect the sensor to the O<sub>2</sub> Sensor-Master.

As soon as the sensor is connected to the O<sub>2</sub> Sensor-Master, the polarization function is automatically activated. Please note: if the sensor was disconnected from the transmitter for longer than 5 minutes, the sensor must be polarized first (polarizing time see section 5.1) to get representative test results.

- **Battery Check:**

Push the switch to the left. If the battery is ok and the O<sub>2</sub> Sensor-Master is operational the green LED lights up. Otherwise, please consult the instruction manual of the O<sub>2</sub> Sensor-Master.

- **Sensor Check:**

For this test the O<sub>2</sub> sensor must be **fully polarized** and the **membrane of the sensor must be dry and clean**.

Expose the sensor (connected to O<sub>2</sub> Sensor-Master) to the air. By pushing the switch to the right to the position «2 – Sensor check», the O<sub>2</sub> Sensor-Master checks whether the electrode current for air measurement delivered by the sensor is within the admissible range (40 to 100 nA for InPro 6800). If the green LED lights up the current for air measurement is within the admissible range.

If the LED does not light up, you should check the battery of the O<sub>2</sub> Sensor-Master (see instruction manual «Accessories»). If the battery is working, than there is probably a problem with your sensor. You should change the electrolyte and/or the membrane body of your sensor. If after a membrane change the LED still does not light up, this means that there is maybe something wrong with the interior body of the sensor. You should than change it (see «Section 6.2»).

**Important!** The sensor check function only verifies the correctness of the electrode current for air measurement. In order to be absolutely sure of the functionality of the sensor, the residual signal in an oxygen free medium should also be controlled (see «Section 6.1.3»).

### 6.1.3 Testing the sensor via a transmitter

A periodic zero current measurement (**no zero point calibration**) is recommended for verification of proper sensor function.

**At the time you carry out the zero current measurement, the sensor must be polarized.**

Zero current measurement can be done by using zeroing gel (order no. 34 100 1032) or nitrogen (N<sub>2</sub>) or carbon dioxide (CO<sub>2</sub>) calibration gases with a purity of at least 99.995 %, alternatively in a sample medium saturated with one of these gases.

After 2 minutes in an oxygen-free sample medium, the reading on the transmitter should drop to below 10 % of the reading in ambient air, and within 10 minutes the value should have dropped to below 1 %.

If the measured values are too high, this suggests a depleted electrolyte or a defective membrane. In the first instance replace the electrolyte, and in the second case exchange both the membrane body and the electrolyte accordingly.

If after such procedures the above mentioned values are still not reached, replace the interior body. If this doesn't solve the problem too send the sensor to your local METTLER TOLEDO representative for inspection.

Many sample media contain volatile substances which, even at very low concentrations, have a clearly perceptible smell. Similarly to oxygen, these sub-stances are able to invade the electrolyte through the gas-permeable membrane. Accordingly, they become noticeable when changing the electrolyte. In most cases, such substances have absolutely no influence on the measuring properties of the sensor. Slight discoloration of the electrolyte also has no effect on the measuring properties.

## 6.2 Changing the electrolyte, the membrane body or the interior body

METTLER TOLEDO DO sensors are supplied with fitted membrane body and have been checked for proper function.

However, if a sensor is to be stored for several months, the electrolyte should be replaced before use.

If the membrane exhibits signs of failure (long response time, increased current in an oxygen-free medium, mechanical damage, etc.) the membrane body has to be replaced.

**Warning!** The O<sub>2</sub> electrolyte has an alkaline pH value of 13. Contact of electrolyte with mucous membrane or eyes is to be avoided. **Therefore protective gloves and safety glasses have to be worn for the following dismantling works.**

If such contact occurs, the affected area should be well rinsed with water. In the case of accident, or should ever any adverse signs appear, get immediate medical attention.

When changing the electrolyte, the membrane body or the interior body, please observe the following instructions (see also the following illustration):

**Attention!** Make sure that this maintenance step is carried out in **clean place**.

1. Unscrew the cap sleeve from the sensor shaft and carefully pull it off the sensor.
2. Pull off the membrane body from the interior body. If it is tight-fitted, eject by pushing it with the flat finger tip. Before electrolyte is refilled, the membrane body must be removed from the cap sleeve.
3. Rinse the interior body with demineralized water and carefully dab it dry with a paper tissue.

**Note: steps 4 to 7 may only be carried out when changing the interior body.**


4. Unscrew the retainer nut of the interior body with an adjustable wrench or with a 3/8" wrench.
5. Remove the interior body by pulling it out of the sensor shaft. If necessary use a plier.

**Warning! Do not twist the interior body.** Otherwise the connection pins can be damaged.

6. Insert the new interior body in the sensor shaft. Turn the interior body in the shaft until the slit of the interior body is aligned with the pin placed in the shaft.
7. Press the body in the shaft and screw the new retainer nut in place.
8. Examine the O-rings visually for mechanical defects, and replace if necessary.
9. Half-fill the new membrane body with O<sub>2</sub> electrolyte and make sure that all bubbles are removed. Air bubbles can be removed by carefully tapping on the membrane body.


10. Slip the membrane body over the interior body while holding the sensor in a vertical position. The excess electrolyte will be displaced and have to be removed with a paper tissue.

**Attention!** No electrolyte, sample media or contamination may be present between the membrane body and the cap sleeve. Please check carefully.

11. Carefully slip the cap sleeve over the fitted membrane body, holding the sensor in a vertical position and screw it tight. The cap sleeve must be clean and dry.
12. After each exchange of electrolyte or membrane body, the sensor has to be repolarized and recalibrated.

## 7 Storage

The sensor can be stored for several months filled with O<sub>2</sub> electrolyte (order no. 34 100 2016) as long as the protective cap is in place. In order to avoid the 6-hour polarization time before reuse of the sensor, it can be kept stored connected to the METTLER TOLEDO O<sub>2</sub> Sensor-Master InPro 6800.

**Attention!** The electrolyte should be replaced before use if the storage period exceeds 3 months.

**If the storage period exceeds 6 months**, the sensor should be **stored dry**, i.e. without any electrolyte in the membrane body.

**Attention!** A sensor being stored dry (without electrolyte in the membrane body) **may on no account be connected to the O<sub>2</sub> Sensor-Master InPro 6800**.

## 8 Product specification

### 8.1 Certificates

Each sensor is delivered with a set of **3.1 B certificates** (complying with EN 10204.3/1.B).

All wetted metal parts (sensor shaft, cap sleeve and membrane body) are identified with a engraved symbol corresponding to the heat number on the paper certificate delivered with the sensor.

Each wetted metal part (sensor shaft, cap sleeve and membrane body) is polished in order to get a surface roughness lower than 0.4 µm (16 µin). This represents a roughness grade number of N5 (according to ISO 1320:1992).

## 8.2 Specifications

**Attention!** Please refer to the supplementary sheet 52 201 146 if the sensor, together with the membrane body T-6800 gas, is been used for measurement of oxygen partial pressure in gases. This supplementary sheet is only included in the packaging of T-6800 gas membrane bodies and kits.


| <b>InPro 6800</b>  |  |
|--|--|
| <b>Measurement principle</b>                               | Amperometric/Polarographic   |
| <b>Working conditions</b>                                  |  |
| Pressure resistance measurement                            | 0.2 – 6 bar absolute<br>[2.9-87.0 psi absolute]                          |
| Mechanical pressure resistance                             | max. 12 bar absolute<br>[174.0 psi absolute]                             |
| Measuring temperature range                                | 0...80 °C [32...176 °F]  |
| Temperature range  | -5...140 °C [23...284 °F]<br>(sterilizable/autoclavable)                 |
| <b>Construction</b>  |  |
| Temperature compensation                                   | Automatic with built-in RTD  |
| Cable connection   | VarioPin (IP68) straight or angled<br>K8S straight (digital ISM sensors) |
| O-ring material  | Silicone FDA and<br>USP VI approved                                      |
| Membrane material  | PTFE/Silicone/PTFE<br>(reinforced with steel mesh)                       |
| Wetted metal parts   | Stainless steel<br>Special material on request                           |
| Surface roughness of<br>wetted metal parts (ISO 1320:1992) | N5 (RA < 0.4 µm [16 µin])  |
| Quick disconnect interior body                             | Standard   |
| Cathode  | Pt   |
| Anode  | Ag   |
| Guard ring   | No   |
| <b>Dimensions</b>  |  |
| Sensor diameter  | 12 or 25 mm (0.47" or 0.98")   |
| Immersion length (α) for 12mm sensor                       | 70, 120, 220, 320, 420 mm<br>(2.8, 4.7, 8.66, 12.6, 16.54")              |
| Immersion length (α) for 25mm sensor                       | 80, 160, 260, 360 mm<br>(3.15, 6.3, 10.24, 14.17")                       |
| <b>Performances</b>  |  |
| Detection limit  | 6 ppb  |
| Accuracy   | " ±[1 % + 6 ppb]<br>of reading in liquids                                |
| Response time at 25°C/77°F (air → N <sub>2</sub> )         | 98 % of final value < 90 s   |
| Sensor signal in ambient air (25 °C/77 °F)                 | 50...110 nA  |
| Residual signal in oxygen-free medium                      | < 0.1 %<br>of the signal in ambient air                                  |
| Maximum flow error   | " 5%   |
| <b>Certification</b>                                       |  |
| EHEDG, 3A  | Yes  |
| 3.1 B (EN 10204.3/1.B)                                     | Yes  |
| ATEX certificate   | Yes  |
| FM Approval  | Yes  |
| FDA/USP VI   | Yes  |
| Quality control  | Yes  |
| <b>Compatibility</b>                                       |  |
| with METTLER TOLEDO transmitters                           | see «Section 9.4»  |
| with METTLER TOLEDO housings                               | see «Section 9.5»  |

## 9 Ordering information

For more detailed information refer to the technical data sheet. Ask your local distributor.

### 9.1 Sensors

#### Intelligent Sensor Management (ISM) O<sub>2</sub> sensor configuration

Example configuration:

InPro 6800/12/120/ISM  
→ Sensor **with** ISM functionality

InPro 6810/25/260  
→ Sensor **without** ISM functionality

#### Sensor without ISM functionality

| Immersion length (a) | 12 mm diameter |            | 25 mm diameter |            |
|----------------------|----------------|------------|----------------|------------|
|                      | straight       | angled     | straight       | angled     |
| 70 mm (2.8")         | 52 200 964     | 52 200 969 | —              | —          |
| 80 mm (3.2")         | —              | —          | 52 200 974     | 52 200 978 |
| 120 mm (4.7")        | 52 200 965     | 52 200 970 | —              | —          |
| 160 mm (6.3")        | —              | —          | 52 200 975     | 52 200 979 |
| 220 mm (8.7")        | 52 200 966     | 52 200 971 | —              | —          |
| 260 mm (10.2")       | —              | —          | 52 200 976     | 52 200 980 |
| 320 mm (12.6")       | 52 200 967     | 52 200 972 | —              | —          |
| 360 mm (14.2")       | —              | —          | 52 200 977     | 52 200 981 |
| 420 mm (16.5")       | 52 200 968     | 52 200 973 | —              | —          |

### 9.2 Accessories

| Accessories                             | Order no.   |
|---|-------------|
| O <sub>2</sub> Sensor-Master InPro 6800 | 52 200 892  |
| O <sub>2</sub> Sensor-Simulator         | 52 200 891  |
| Adapter T82 socket – VP plug            | 52 200 939  |
| Adapter VP socket – T82 plug            | 52 200 940  |
| VP cable VP6-ST/3 m                     | 52 300 108  |
| VP cable VP6-ST/5 m                     | 52 300 109  |
| VP cable VP6-HT/3 m                     | 52 300 112  |
| VP cable VP6-HT/5 m                     | 52 300 113  |
| VP cable VP8-ST/3 m                     | 52 300 354  |
| VP cable VP8-ST/5 m                     | 52 300 355  |
| VP cable VP8-HT/3 m                     | 52 300 361  |
| VP cable VP8-HT/5 m                     | 52 300 362  |
| Zeroing gel                             | 34 100 1032 |

For other cable lengths or types, please contact your local Mettler-Toledo representative.

### 9.3 Spare parts

| Spare parts   | Order no.   |
|---|-------------|
| Membrane body, single T-96  | 52 200 071  |
| Membrane kit T-96<br>(4 membrane bodies, 1 O-ring set,<br>25 ml Electrolyte bottle) | 52 200 024  |
| Membrane bodies (20 pcs.) T-96  | 52 200 791  |
| Electrolyte bottle (25 ml)  | 34 100 2016 |
| Interior body InPro 6800<br>(with quick disconnect)                                 | 52 200 899  |
| Cap sleeve N-Type   | 52 200 037  |
| Cap sleeve with protective cage P-Type  | 52 200 038  |

### 9.4 Recommended transmitters

| Transmitter                               | Order no.  |
|---|--|
| O <sub>2</sub> Transmitter M 700          | Ask your local distributor<br>or refer to the technical<br>data sheet. |
| O <sub>2</sub> Transmitter 4100 e FF      |  |
| O <sub>2</sub> Transmitter 4100 e/2 (x) H |  |
| O <sub>2</sub> Transmitter 4100 PA        |  |
| O <sub>2</sub> Transmitter 4100 e         |  |
| O <sub>2</sub> Transmitter 4050e          |  |

### 9.5 Recommended housings

| Housings                          | 12 mm Ø | 25 mm Ø |   |
|-----------------------------------|---------|---------|---|
| <b>Static housings</b>            |         |         | Please ask<br>your local<br>Mettler-Toledo<br>representative. |
| INGOLD "safety<br>weld-in socket" | –       | ✓       |   |
| InFit 761 CIP                     | ✓       | –       |   |
| <b>Retractable housings</b>       |         |         |   |
| InFit 777 e                       | ✓       | –       |   |
| InFit 797 e                       | ✓       | –       |   |
| <b>Immersion housing</b>          |         |         |   |
| InDip 550                         | ✓       | –       |   |

**Note:** The housings are available in different versions. Please contact your distributor to get the right ordering information.

## 10 Theory of the polarographic sensor

### 10.1 Introduction

Two types of electrodes are employed in analytical work: **potentiometric** and **amperometric** electrodes.


- Potentiometric electrodes develop a voltage generated by the activity of a particular ion. Examples of such electrodes are glass electrodes (like pH electrodes) and most ion-selective electrodes. Their individual potentials cannot be determined. The measurable quantity is the difference of potential between the measuring electrode and an inert reference electrode. The potential of the reference electrode must be constant.

All potentiometric electrodes are subject to Nernst's law and for this reason electrodes and measuring instruments are in most cases interchangeable. An important requirement of potentiometric measurements is the virtually currentless determination of the electrode voltage. During measurement no chemical reaction occurs and the solution remains in equilibrium.


- In the case of **amperometric electrodes**, such as the **oxygen electrode**, activity measurement is based on a current measurement.

**The oxygen electrode consists of a cathode and an anode conductively connected by an electrolyte.** A suitable polarization voltage between the anode and the cathode selectively reduces the oxygen at the cathode.

Cathode reaction:

Anode reaction:

These chemical reactions result in an electric current which is proportional to the oxygen partial pressure (PO<sub>2</sub>). The oxygen electrode consumes oxygen which is continuously extracted from the solution. The viscosity and flow rate of the solution are therefore important parameters.

The electrode current of an oxygen electrode is determined not only by the oxygen partial pressure but by many other electrode parameters. The electrode currents of different electrode types may differ by several powers of ten. For this reason oxygen electrodes and amplifiers cannot be freely interchanged.


## 10.2 Principle of the design of an oxygen electrode

There are two main types of oxygen electrodes:

- Electrodes **without** membrane
- Electrodes **with** gas-permeable membrane (Clark Principle)

The membrane electrode according to Clark is today most widely used. As compared to the electrode without membrane it possesses the following advantages:

- Oxygen measurement in gases and solutions
- No mutual contamination of electrode and solution
- No or little dependency on flow

In the case of the Clark electrode, geometrical configuration is very important. In particular, the thickness of the electrolyte film between the cathode and the membrane must be within narrow tolerances so as to ensure good linearity and a low zero current (current in a nitrogen atmosphere). The figure above exemplifies the principle of the design of a Clark-type oxygen electrode.

## 10.3 Parameters determining current

The quantity of oxygen diffused in and the magnitude of the electrode current are influenced by the following parameters:

- Oxygen partial pressure of the solution
- Membrane material and thickness
- Size of cathode
- Polarization voltage
- Temperature
- Flow conditions in the solution


Fick's law gives the mathematical relationship between these parameters:

$$I = k \cdot D \cdot a \cdot A \cdot \frac{pO_2}{X}$$

- I = electrode current  
 k = constant  
 D = coefficient of diffusion of O<sub>2</sub> in the membrane  
 a = solubility of O<sub>2</sub> in the membrane material  
 A = cathode surface  
 pO<sub>2</sub> = oxygen partial pressure in the solution  
 X = thickness of gas-permeable membrane

## 10.4 Polarization voltage

The voltage between the anode and the cathode is so selected that the oxygen is fully (> A, see polarogram) reduced while other gases are unaffected (< D). The ideal voltage for the Pt/Ag/AgCl system is between -500 and -750 mV.

The polarization voltage should remain as constant as possible. Besides a constant source of voltage the following prerequisites must be met:

The electrical resistance of electrolyte film must not exceed a particular value so as to avoid a voltage drop. The anode must have a large surface so as to prevent polarization of the anode by the electrode current.

## 10.5 Temperature

The temperature dependence of the current passing through an oxygen electrode when referred to a constant O<sub>2</sub> partial pressure is determined mainly by the properties of the gas-permeable membrane.

## 10.6 Dependence on flow

With most oxygen electrodes the electrode current is smaller in stagnant solutions than in agitated ones. The consumption of oxygen by the electrode results in extraction of oxygen from the solution in the close vicinity of the cathode outside the membrane. The oxygen is replaced by diffusion. If the electrode current is strong, the solution cannot fully restore the oxygen by diffusion. This results in an electrode current weaker than would correspond to conditions in the solution. In agitated solutions the oxygen is conveyed to the surface of the membrane not only by diffusion but additionally by the flow (convection). In that case no oxygen impoverishment occurs at the membrane surface.

A high degree of flow dependence occurs mainly with large cathodes, thin and highly permeable membranes, i.e. where electrode currents are large. The problem of flow dependence is often solved by prescribing a minimum flow rate.

In METTLER TOLEDO InPro 6800 electrodes, the thin PTFE membrane determining the electrode current (i.e. the actual measuring signal) is separated from the sample solution by a relatively thick silicone membrane. This latter is highly permeable to oxygen molecules and thus acts as an oxygen reservoir. The diffusion of oxygen out of the sample solution into the silicone membrane is spread over a wide area. Since this results in less oxygen being extracted from the sample solution per unit area, the double PTFE/Silicone membrane forms an effective buffer against disturbances due to hydrodynamic flow.

## 10.7 Oxygen partial pressure – oxygen concentration

The electrode current depends on the oxygen partial pressure and the O<sub>2</sub> permeability of the membrane – but not on the O<sub>2</sub> solubility in the solutions. The oxygen concentration in mg O<sub>2</sub>/l (CL) cannot therefore be determined directly with an electrode.

According to Henry's law the oxygen concentration is proportional to its partial pressure (PO<sub>2</sub>).

$$CL = pO_2 \cdot a$$

$a$  = solubility factor

If « $a$ » is constant, the oxygen concentration can be determined by means of the electrode. This applies at constant temperature and with dilute aqueous solutions such as drinking-water.

The solubility factor is strongly influenced not only by the temperature but also by the composition of the solution:

| Medium, sat. with air | Solubility at 20°C (68°F) and 760 mm Hg |
|-----------------------|---|
| Water                 | 9.2 mg O <sub>2</sub> /l                |
| 4 mol/l KCl           | 2 mg O <sub>2</sub> /l                  |
| 50 % Methanol-water   | 21.9 mg O <sub>2</sub> /l               |

Although the solubilities vary widely, the oxygen electrode gives the same reading in all three solutions. Thus, determination of the oxygen concentration is only possible with constant and known solubility factors « $a$ ».

Solubility may be determined by a Winkler titration or the method developed by Käppeli and Fiechter.

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- CH** **Mettler-Toledo (Schweiz) GmbH**, Im Langacher, Postfach, CH - 8606 Greifensee, Switzerland
- D** **Mettler-Toledo GmbH**, Prozeßanalytik, Ockerweg 3, D - 35396 Gießen, Germany
- F** **Mettler-Toledo Analyse Industrielle S.A.S.**, 30 Bld. de Douaumont, F - 75017 Paris, France,
- USA** **Mettler-Toledo Ingold, Inc.**, 36 Middlesex Turnpike, Bedford, MA 01730, USA